

PHLEBOTOMY ORDER OF DRAW							
FIRST							LAST
BLOOD CULTURE BOTTLES Must collect with syringe. Use sterile technique for site prep.	CITRATE	SST	PLAIN RED	HEPARIN	EDTA	FLUORIDE	ACD
	EXCYTE ESR		NAVY RED LABEL		NAVY LAV LABEL		
	Observe fill lines on light blue top tubes						

CAPILLARY ORDER OF DRAW					
FIRST					LAST
CBG'S	EDTA	HEPARIN	PLAIN RED	SST	NEONATAL